

What Does It Measure?

The Visioline® VL 650 is the ideal instrument to objectively analyze the deeper lines and **macro wrinkles** such as crow's feet. It is a further development of the renowned Quantirides® system.

The Measuring Principle





The measurement is based on **skin replica and oblique lighting**. The replica is illuminated at an angle of 35° and the mountains representing the wrinkles of the skin produce measurable shadows. They are digitalized by a **high resolution** camera mounted vertically to the replica and serve as a basis for different arithmetical calculations (length, depth and area of the wrinkles in µm).

Fields of Application

- This economical system is indispensable in **efficacy testing and claim support** for cosmetic anti-wrinkle products.
- Perfect tool for **multicentric studies**, as the replicas can be collected over a long time in different places and then be evaluated together.

Advantages

- **Easy to make replicas** in all sizes from all different body sites.
- On the **mount** the replica can be moved very accurately in x and y direction by screws.
- To analyse the same site before and after treatment, perfect placing is possible with the help of a histogram of the shadows and an **overlay mask** of the previous replica.
- Measurements and calibration can be performed **very quickly**.
- **3D and false colours** for impressive marketing purposes available.
- All study data are automatically stored in a **database** in the software.



Technical Data

Power supply: illumination: external, 12 VDC, Camera: USB; Dimensions: 15.3 x 21.7 x 21.7 cm;
 Measurement area: from 13.5 x 18 mm to 16.5 x 22 mm; Weight: 4 kg; Port: USB; Light source: white LED under 35°
 (± 0.5°); Shadow length determination in µm; xy-Stage: Resolution: 1µm, Accuracy: 2µm, Range: 10 mm;
 Camera: 2560 x 1920 Pixel, 5 MPix; Objective: Focal length: 25 mm, Aperture: 1.4 - 16;
 Computer: Windows® 10, performance must meet system requirements, USB 2.0
 Measurement principle: optical, reflected light from replica

Technical changes may be made without prior notice.

Courage+Khazaka electronic GmbH since 1986
 Mathias-Brüggen-Str. 91 · 50829 Köln · GERMANY
 phone +49 221 95 64 99 0 · fax +49 221 95 64 99 1
 info@courage-khazaka.de · www.courage-khazaka.de

