### What Does It Measure?

Easy and quick method to analyse dandruff in number and size.

# **The Principle**


The system consists of a device in which a **petri dish** with the collected dandruff is inserted. A circularly arranged **LED light source** illuminates the sample homogeneously on a dark background. The high resolution camera above the sample takes the image and the **software detects** all dandruff and categorizes it in 9 different size classes.

## **Fields of Application**

Hair care is a large field in the cosmetic industry. Analysing dandruff is one of the **most important concerns** in research of hair care products and their **efficacy testing**.

# **Advantages**

- Quick and easy handling of the device.
- Simple connection to the computer by **USB**
- Dandruff is evaluated by number & size (in pixel and mm²) for 9 different size categories which can be determined by the user.
- The average of up to 4 images is automatically calculated.
- Software conveniently allows evaluation of complete studies.
- Easy **calibration** of the system.





#### Technical Data

Dimensions:  $13.5 \times 13 \times 15.3 \text{ cm}$  (H x W x D), Bevel:  $10^{\circ}$  on the front,  $60^{\circ}$  on the back, Opening:  $9.2 \times 3.1 \text{ cm}$  (W x H), Weight: 1.35 kg, Port: USB 2.0, Power supply: Input: 110-240 V, 50-60 Hz, Output: DC 12V/max. 4A, Light source: white LED light, arranged circularly, USB-Camera: 1/2" CMOS, Resolution:  $1280 \times 1024 \text{ Pixel} = 1.3 \text{ MPixel}$ , max. 25 images/second, Objective: M12; 6 mm focal length, distance camera to sample: approx. 9 cm, Petri dish:  $\emptyset 8.5 \text{ cm}$ , visible field  $\emptyset 7.5 \text{ cm}$ .

Courage+Khazaka electronic GmbH since 1986 Mathias-Brüggen-Str. 91 · 50829 Köln · GERMANY

phone +49 221 95 64 99 0 · fax +49 221 95 64 99 1 info@courage-khazaka.de · www.courage-khazaka.de

